



The Biotechnology Education Company ®

Mapping of Restriction Sites on Plasmid DNA

EDVO-Kit
105

See Page 3 for storage instructions.

EXPERIMENT OBJECTIVE:

The objective of this experiment module is to develop an understanding of the principles of DNA mapping using various restriction enzymes to generate DNA fragments.

EDVOTEK, Inc. • 1-800-EDVOTEK • www.edvotek.com

Mapping of Restriction Sites on Plasmid DNA

Table of Contents

	Page
Experiment Components	3
Experiment Requirements	3
Background Information	4
Experiment Procedures	
Experiment Overview and General Instructions	7
Agarose Gel Electrophoresis	9
Size Determination of DNA Restriction Fragments	10
Mapping of DNA Restriction Sites	12
Study Questions	13
Instructor's Guidelines	
Notes to the Instructor and Pre-Lab Preparations	15
Experiment Results and Analysis	21
Study Questions and Answers	22
Appendices	23
Material Safety Data Sheets	34

All components are intended for educational research only. They are not to be used for diagnostic or drug purposes, nor administered to or consumed by humans or animals.

THIS EXPERIMENT DOES NOT CONTAIN HUMAN DNA. None of the experiment components are derived from human sources.

EDVOTEK, The Biotechnology Education Company, and InstaStain are registered trademarks of EDVOTEK, Inc.. Ready-to-Load, UltraSpec-Agarose and FlashBlue are trademarks of EDVOTEK, Inc.



Experiment Components

DNA samples are stable at room temperature. However, if the experiment will not be conducted within one month of receipt, it is recommended that the DNA samples be stored in the refrigerator.

DNA samples do not require heating prior to gel loading.

READY-TO-LOAD™ DNA SAMPLES FOR ELECTROPHORESIS

- A Standard DNA Fragments
- B Plasmid cut with Enzyme 1
- C Plasmid cut with Enzyme 2
- D Plasmid cut with Enzyme 1 and Enzyme 2

REAGENTS & SUPPLIES

- UltraSpec-Agarose™ powder
- Concentrated electrophoresis buffer
- FlashBlue™ DNA Stain
- InstaStain® Blue cards
- Practice Gel Loading Solution
- 1 ml pipet
- Microtipped Transfer Pipets

Note: If you ordered Experiment #105-Q, the experiment components include InstaStain® Ethidium bromide instead of FlashBlue™ and InstaStain® Blue DNA stains.

Requirements



- Horizontal gel electrophoresis apparatus
- D.C. power supply
- Automatic micropipets with tips
- Balance
- Microwave, hot plate or burner
- Pipet pump
- 250 ml flasks or beakers
- Hot gloves
- Safety goggles and disposable laboratory gloves
- Small plastic trays or large weigh boats (for gel destaining)
- DNA visualization system (white light)
- Distilled or deionized water

Background Information

The Human Genome and other genome projects are extremely significant accomplishments with important applications to biology and medicine. The explosion of this new information is leading to dramatic changes in the way we are able to improve life. Part of the challenge in dealing with the enormous amounts of data is to determine what genes are responsible for different functions. Scientists must determine the location of genes through DNA mapping, and then begin the arduous task of determining what the individual genes do.

Mapping the positions of restriction enzyme cleavage sites on a DNA molecule is an important prerequisite to DNA sequencing, which provides the primary nucleotide sequence information in DNA. Mapping involves the determination of the relative distances between restriction enzyme cleavage sites. An illustrative analogy would be somewhat similar to the following: If DNA mapping were compared to identifying the streets on a city map, then DNA sequencing would be analogous to identifying the specific houses on the streets.

DNA mapping is performed by determining the size of the DNA fragments generated by single or combinations of restriction enzyme digestions, and subsequent construction of a DNA map. For example:



- Consider a 5000 base pair, circular plasmid DNA containing single recognition sites for enzymes A, B, and C (Figure 1).
- Going in a clockwise direction from A, the distances between

A and B is 500
B and C is 1500
C and A is 3000

These assignments are made based on the size of the entire circular plasmid, which is 5,000 base pairs (Figure 2).

To obtain a reference point, the cleavage site at A will be arbitrarily assigned as position zero. All three enzymes will cleave the plasmid once to produce a linear molecule of 5000 base pairs. Different combinations of these enzymes will produce the following DNA fragments (in base pairs):

A+B	A+C	B+C	A+B+C
4500	3000	3500	3000
500	2000	1500	1500
			500



Background Information

This data shows that the cleavage site at B is closest to A since cleavage A+B generated the smallest fragment (500) out of all the pairs of enzymes. The shortest distance between A and C is 2000 base pairs since the smallest fragment in the A + C pair is 2000. Similarly, the shortest distance between B and C is 1500 base pairs. It remains to be determined if B is in between A and C (Figure 1) or alternatively, B is between C and A (going in a clockwise direction from A around the plasmid, Figure 3).

If C was in between A and B, the 500 base pair fragment would have been cleaved into two smaller fragments. However, when all three enzymes are used, the 500 base pair fragment remains. In addition, only the 2000 base pair fragment found in the A + C pair is cleaved into 1500 and 500 base pair fragments when all three enzymes are used, verifying the location of B. This kind of logic enables the construction of a map, as previously shown, from DNA fragment sizes.

Note that the data from this experiment cannot tell us the absolute orientation of the cleavage sites since it can lead to an alternative map as shown in Figure 2. However, the relative positions are still the same (B is in between A and C). The assignment (Figure 1 or Figure 3) can be made upon further analysis.

Unknown DNA fragment sizes are determined by comparing the relative mobilities of DNA fragments of known size as standards. DNA fragments, from plasmid digests, and standard DNA fragments (also known as markers) are electrophoretically separated in parallel on the same agarose gel. After electrophoretic separation, DNA fragments are stained for visualization, and migration distances of known and unknown fragments are measured.




Figure 3:
Construction of a
plasmid map using three
enzymes.

Background Information

Standard fragments are used to make a standard curve by plotting their size on the y-axis versus the migration distance on the x-axis. The size of the fragments on the y-axis are expressed as the log of the number of base pairs they contain or the log of their molecular weight. Most of the plotted data obtained from the markers will yield a straight line. The migration distance of the unknown DNA fragment(s) are located on the X-axis and their size is estimated from the standard curve.


After determining the size of the DNA fragments generated by single and combinations of restriction enzymes, a DNA map is constructed as previously described.

In this experiment, you will determine the relative locations of three restriction enzyme cleavage sites on a circular plasmid DNA. The plasmid has been cleaved with three restriction enzymes. Enzyme 1 cleaves the plasmid once at site A. Assume that the Enzyme 1 site is at position 0. Enzyme 2 and 3 also cut the plasmid once at sites B and C. The objective is to calculate the distances in base pairs between the points of cleavage and to determine whether the Enzyme 1 site is in between the Enzyme 2 sites.

Quick Reference:

A standard curve will be made on semi-log graph paper. The following are the Standard DNA fragment sizes, which are expressed in base pairs.

23130	9416	6557
4361	3000	2322
2027	725	570



Experiment Overview and General Instructions

EXPERIMENT OBJECTIVE:

The objective of this experiment module is to develop an understanding of the principles of DNA mapping using various restriction enzymes to generate DNA fragments.

LABORATORY SAFETY

1. Gloves and goggles should be worn routinely as good laboratory practice.
2. Exercise extreme caution when working with equipment that is used in conjunction with the heating and/or melting of reagents.
3. DO NOT MOUTH PIPET REAGENTS - USE PIPET PUMPS.
4. Exercise caution when using any electrical equipment in the laboratory.
5. Always wash hands thoroughly with soap and water after handling reagents or biological materials in the laboratory.



LABORATORY NOTEBOOK RECORDINGS:

Address and record the following in your laboratory notebook or on a separate worksheet.

Before starting the Experiment:

- Write a hypothesis that reflects the experiment.
- Predict experimental outcomes.

During the Experiment:


- Record (draw) your observations, or photograph the results.

Following the Experiment:

- Formulate an explanation from the results.
- Determine what could be changed in the experiment if the experiment were repeated.
- Write a hypothesis that would reflect this change.

Experiment Overview: Flow Chart

Experiment Procedure



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 1989, 1992, 1994, 1997, 1998, 2000, 2004, 2007, 2009, EDVOTEK, Inc., all rights reserved. EVT 100202AM

Agarose Gel Electrophoresis

Prepare the Gel

1. Prepare an agarose gel with specifications summarized below. Your instructor will specify which DNA stain you will be using.



- Agarose gel concentration required: 0.8%
- Recommended gel size: 7 x 7 cm or 7 x 14 cm (two gels)
- Number of sample wells required: 4
- Placement of well-former template: first set of notches (7 x 7 cm)
first & third set of notches (7 x 14 cm)

For gels to be stained with FlashBlue™ or InstaStain® Blue, prepare gels according to Appendix A.

For gels to be stained with InstaStain® Ethidium bromide, prepare gels according to Appendix B.

Step-by-step guidelines for agarose gel preparation are summarized in Appendix D.

Load the Samples

2. Load the DNA samples in tubes A - D into the wells in consecutive order.
 - For gels to be stained with FlashBlue™ or InstaStain® Blue, fill wells with 35 - 38 µl.
 - For gels to be stained with InstaStain® Ethidium Bromide, fill wells with 18 - 20 µl.

Lane	Tube	Sample
1	A	Standard DNA Fragments
2	B	Plasmid cut with Enzyme 1
3	C	Plasmid cut with Enzyme 2
4	D	Plasmid cut with Enzyme 1 and Enzyme 2


Run the Gel

3. After DNA samples are loaded, connect the apparatus to the D.C. power source and set the power source at the required voltage.
4. Check that current is flowing properly - you should see bubbles forming on the two platinum electrodes. Conduct electrophoresis for the length of time specified by your instructor.
5. After electrophoresis is completed, proceed to DNA staining and visualization. Refer to Appendix E, F, G, or H for the appropriate staining instructions.
6. Document the results of the gel by photodocumentation.

Alternatively, place transparency film on the gel and trace it with a permanent marking pen. Remember to include the outline of the gel and the sample wells in addition to the migration pattern of the DNA bands.

Reminders:


During electrophoresis, the DNA samples migrate through the agarose gel towards the positive electrode. Before loading the samples, make sure the gel is properly oriented in the apparatus chamber.



This exercise focuses on the first step for mapping DNA restriction sites, which is to determine the size of "unknown" DNA fragments generated after electrophoresis. The assignment of sizes for DNA fragments separated by agarose gel electrophoresis can have $\pm 10\%$ margin of error. The sizes of the "unknowns" will be extrapolated by their migration distances relative to the Standard DNA Fragments (Sample A), for which the fragment sizes are known.

1. Measure and record the distance traveled in the agarose gel by each Standard DNA fragment (except the largest 23,130 bp fragment, which will not fit in a straight line in step 4).

Figure 1



In each case, measure from the lower edge of the sample well to the lower end of each band. Record the distance traveled in centimeters (to the nearest millimeter).

2. Label the semi-log graph paper:
 - A. Label the non-logarithmic horizontal x-axis "Migration Distance" in centimeters at equal intervals.
 - B. Label the logarithmic vertical y-axis "Log base pairs". Choose your scales so that the data points are well spread out. Assume the first cycle on the y-axis represents 100-1,000 base pairs and the second cycle represents 1,000-10,000 base pairs.
3. For each Standard DNA fragment, plot the measured migration distance on the x-axis versus its size in base pairs, on the y-axis.
4. Draw the best average straight line through all the points. The line should have approximately equal numbers of points scattered on each side of the line. Some points may be right on the line (see Figure 1 for an example).
5. Measure the migration distance of each of the "unknown" fragments from samples B, C, and D.


Quick Reference:

Standard DNA fragment sizes - length is expressed in base pairs.

23130	9416	6557
4361	3000	2322
2027	725	570

6. Using the graph of the Standard DNA fragments, determine the sizes in base pairs of each "unknown" fragment.
 - Find the migration distance of the unknown fragment on the x-axis. Draw a vertical line from that point until the standard graph line is intersected.
 - From the point of intersection, draw a second line horizontally to the y-axis and determine the approximate size of the fragment in base pairs (refer to Figure 1 for an example).







Mapping of DNA Restriction Sites

The size of the plasmid used in this experiment is 4300 bp.


1. Draw a circle representing a 4300 bp plasmid on a transparent sheet of acetate.
2. Mark the positions of Enzyme #2 (Lane 3) sites corresponding to the sizes of fragments obtained upon digestion of the plasmid on the gel.
3. Draw a second circle representing a 4300 bp plasmid on a transparent sheet of acetate.
4. Mark the position of the Enzyme #1 (Lane 2) site at the top (12:00 o'clock).




5. To draw a composite map of both enzymes, overlay the Enzyme #2 map on top of the Enzyme #1 map.
6. Keeping the Enzyme #1 site at the 12:00 o'clock position, rotate the Enzyme #2 map until the relative distances between the sites approximate the relative sizes of the fragments of Enzyme #1 and #2 combined.
7. Specify, in base pairs, the distances between all the sites.

 Material Safety Data Sheet May be used to comply with OSHA's Hazard Communication Standard, 29 CFR 1910.1200. Standard must be consulted for specific requirements.	
IDENTITY (As Used on Label and List) Agrose	
Section I - Manufacturer's Name EDVOTEK, Inc. Address (Number, Street, City, State, Zip Code) 14676 Rothgeb Drive Rockville, MD 20850	
Section II - Hazardous Ingredients/Identify Information Hazardous Components Specific Chemical Identity, Common Name(s) OSHA PEL ACGIH TLV Recommended % (Optional) Standard This product contains no hazardous materials as defined by the OSHA Hazard Communication Standard.	
Section III - Physical/Chemical Characteristics Boiling Point For 1% solution (94 F) Specific Gravity (H ₂ O = 1) No data Vapor Pressure (mm Hg) No data Melting Point No data Vapor Density (AIR = 1) No data Evaporation Rate (Butyl Acetate = 1) No data Solubility in Water Insoluble - cold	
Appearance and Odor White powder, no odor	
Section IV - Physical/Chemical Characteristics N.D. = No data Flash Point (Method Used) No data Flammable Limits LEL N.D. UEL N.D. Extinguishing Media Water spray, dry chemical, carbon dioxide, halon or standard foam	
Special Fire Fighting Procedures Possible fire hazard when exposed to heat or flame	
Unusual Fire and Explosion Hazards None	

Section V - Reactivity Data Stability Unstable Stable X None Incompatibility No data available Hazardous Decomposition or Byproducts	
Hazardous Polymerization May Occur Will Not Occur X None Section VI - Health Hazard Data Routes of Entry: Inhalation? Yes Skin? Yes Ingestion? Yes Health Hazards (Acute and Chronic) Ingestion: Large amounts may cause diarrhea Carcinogenicity: No data available IARC Monographs? OSHA Regulation?	
Signs and Symptoms of Exposure No data available Medical Conditions Generally Aggravated by Exposure No data available Emergency First Aid Procedures Treat symptomatically and supportively	
Section VII - Precautions for Safe Handling and Use Steps to be Taken in case of Material Released or Spilled Sweep up and place in suitable container for disposal Waste Disposal Method Normal solid waste disposal Precautions to be Taken in Handling and Storing None	
Other Precautions None	
Section VIII - Control Measures Respiratory Protection (Specify Type) Chemical cartridge respirator with full facepiece. Ventilation Local Exhaust Yes Special Mechanical (General) Yes Other None Protective Gloves Yes Eye Protection Splash proof goggles Other Protective Clothing or Equipment Impervious clothing to prevent skin contact Work/Hygiene Practices None	

 Material Safety Data Sheet May be used to comply with OSHA's Hazard Communication Standard, 29 CFR 1910.1200. Standard must be consulted for specific requirements.	
IDENTITY (As Used on Label and List) 50x Electrophoresis Buffer	
Section I - Manufacturer's Name EDVOTEK, Inc. Address (Number, Street, City, State, Zip Code) 14676 Rothgeb Drive Rockville, MD 20850	
Section II - Hazardous Ingredients/Identify Information Hazardous Components Specific Chemical Identity, Common Name(s) OSHA PEL ACGIH TLV Recommended % (Optional) Standard This product contains no hazardous materials as defined by the OSHA Hazard Communication Standard.	
Section III - Physical/Chemical Characteristics Boiling Point No data Specific Gravity (H ₂ O = 1) No data Vapor Pressure (mm Hg) No data Melting Point No data Vapor Density (AIR = 1) No data Evaporation Rate (Butyl Acetate = 1) No data Solubility in Water Appreciable, (greater than 10%)	
Appearance and Odor Clear, liquid, slight vinegar odor	
Section IV - Physical/Chemical Characteristics N.D. = No data Flash Point (Method Used) No data Flammable Limits LEL N.D. UEL N.D. Extinguishing Media Use extinguishing media appropriate for surrounding fire.	
Special Fire Fighting Procedures Wear protective equipment and SCBA with full facepiece operated in positive pressure mode.	
Unusual Fire and Explosion Hazards None identified	

Section V - Reactivity Data Stability Unstable Stable X None Incompatibility Strong oxidizing agents Hazardous Decomposition or Byproducts Carbon monoxide, Carbon dioxide	
Hazardous Polymerization May Occur Will Not Occur X None Section VI - Health Hazard Data Routes of Entry: Inhalation? Yes Skin? Yes Ingestion? Yes Health Hazards (Acute and Chronic) None Carcinogenicity: None identified IARC Monographs? OSHA Regulation?	
Signs and Symptoms of Exposure Irritation to upper respiratory tract, skin, eyes Medical Conditions Generally Aggravated by Exposure None Emergency First Aid Procedures Ingestion: If conscious, give large amounts of water Eyes: Flush with water. Inhalation: Move to fresh air. Skin: Wash with soap and water	
Section VII - Precautions for Safe Handling and Use Steps to be Taken in case of Material Released or Spilled Mop up spill and rinse with water, or collect in absorbent material and dispose of the absorbent material Waste Disposal Method Dispose in accordance with all applicable federal, state, and local environmental regulations. Precautions to be Taken in Handling and Storing Avoid eye and skin contact.	
Other Precautions None	
Section VIII - Control Measures Respiratory Protection (Specify Type) None Ventilation Local Exhaust Yes Special Mechanical (General) Yes Other None Protective Gloves Yes Eye Protection Safety goggles Other Protective Clothing or Equipment None Work/Hygiene Practices None	


 Material Safety Data Sheet May be used to comply with OSHA's Hazard Communication Standard, 29 CFR 1910.1200. Standard must be consulted for specific requirements.	
IDENTITY (As Used on Label and List) Practice Gel Loading Solution	
Section I - Manufacturer's Name EDVOTEK, Inc. Address (Number, Street, City, State, Zip Code) 14676 Rothgeb Drive Rockville, MD 20850	
Section II - Hazardous Ingredients/Identify Information Hazardous Components Specific Chemical Identity, Common Name(s) OSHA PEL ACGIH TLV Recommended % (Optional) Standard This product contains no hazardous materials as defined by the OSHA Hazard Communication Standard.	
Section III - Physical/Chemical Characteristics Boiling Point No data Specific Gravity (H ₂ O = 1) No data Vapor Pressure (mm Hg.) No data Melting Point No data Vapor Density (AIR = 1) No data Evaporation Rate (Butyl Acetate = 1) No data Solubility in Water Soluble	
Appearance and Odor Blue liquid, no odor	
Section IV - Physical/Chemical Characteristics Flash Point (Method Used) No data Flammable Limits LEL No data UEL No data Extinguishing Media Dry chemical, carbon dioxide, water spray or foam	
Special Fire Fighting Procedures Use agents suitable for type of surrounding fire. Keep upwind, avoid breathing hazardous sulfur oxides and bromides. Wear SCBA.	
Unusual Fire and Explosion Hazards Unknown	

Section V - Reactivity Data Stability Unstable Stable X None Incompatibility None Hazardous Decomposition or Byproducts Sulfur oxides, and bromides	
Hazardous Polymerization May Occur Will Not Occur X None Section VI - Health Hazard Data Routes of Entry: Inhalation? Yes Skin? Yes Ingestion? Yes Health Hazards (Acute and Chronic) Acute eye contact: May cause irritation. Carcinogenicity: No data available IARC Monographs? OSHA Regulation?	
Signs and Symptoms of Exposure May cause skin or eye irritation Medical Conditions Generally Aggravated by Exposure None reported Emergency First Aid Procedures Treat symptomatically and supportively. Rinse contacted area with copious amounts of water.	
Section VII - Precautions for Safe Handling and Use Steps to be Taken in case of Material Released or Spilled Wipe eye and skin protection and mop spill area. Rinse with water. Waste Disposal Method Observe all federal, state, and local regulations. Precautions to be Taken in Handling and Storing Avoid eye and skin contact.	
Other Precautions None	
Section VIII - Control Measures Respiratory Protection (Specify Type) None Ventilation Local Exhaust Yes Special Mechanical (General) Yes Other None Protective Gloves Yes Eye Protection Splash proof goggles Other Protective Clothing or Equipment None required Work/Hygiene Practices Avoid eye and skin contact	


Material Safety Data Sheets

Full-size (8.5 x 11") pdf copy of MSDS is available at www.edvotek.com or by request.

105
Experiment

 Material Safety Data Sheet May be used to comply with OSHA's Hazard Communication Standard, 29 CFR 1910.1200. Standard must be consulted for specific requirements.			
IDENTITY (As Used on Label and List) InstaStain® Ethidium Bromide			
<small>Note: Blank space for information that is not applicable, or no information is available, this space must be marked to indicate that.</small>			
Section I Manufacturer's Name InstaStain, Inc. P.O. Box 1232 West Bethesda, MD 20827			
Emergency Telephone Number (301) 251-5990 Telephone Number for information (301) 251-5990 Date Prepared 10/05/06 Signature of Preparer (optional)			
Section II - Hazardous Ingredients/Identify Information Hazardous Components Specific Chemical Identity, Common Name(s) OSHA PEL ACGIH TLV Other Limits Recommended % (Optional)			
Ethidium Bromide Data not available (2,7-Diamino-10-Ethyl-9-Phenylphenanthridinium Bromide)			
CAS # 139-33-3			
Section III - Physical/Chemical Characteristics			
Boiling Point	No data	Specific Gravity (H ₂ O = 1)	No data
Vapor Pressure (mm Hg.)	No data	Melting Point	No data
Vapor Density (AIR = 1)	No data	Evaporation Rate (Butyl Acetate = 1)	No data
Solubility in Water Soluble			
Appearance and Odor Chemical bound to paper, no odor			
Section IV - Physical/Chemical Characteristics			
Flash Point (Method Used)	No data	Flammable Limits	LEL N.D. UEL N.D.
Extinguishing Media Water spray, carbon dioxide, dry chemical powder, alcohol or polymer foam			
Special Fire Fighting Procedures Wear protective clothing and SCBA to prevent contact with skin & eyes			
Unusual Fire and Explosion Hazards Emits toxic fumes			

Section V - Reactivity Data			
Stability	Unstable Stable	Conditions to Avoid	None
Incompatibility	Strong oxidizing agents		
Hazardous Decomposition or Byproducts Carbon dioxide, nitrogen oxide, hydrogen bromide gas			
Hazardous Polymerization	May Occur Will Not Occur	Conditions to Avoid	None
Section VI - Health Hazard Data			
Route(s) of Entry:	Inhalation? Yes	Skin? No	Ingestion? Yes
Health Hazards (Acute and Chronic) May cause eye irritation			
Skin: May cause skin irritation			
Eyes: May cause eye irritation			
Inhalation: Cyanosis			
NTP? IARC Monographs? OSHA Regulation? Meets criteria for proposed OSHA medical records rule PEREAC 47.30420.82			
Signs and Symptoms of Exposure No data available			
Medical Conditions Generally Aggravated by Exposure No data available			
Emergency First Aid Procedures Treat symptomatically and supportively			
Section VII - Precautions for Safe Handling and Use			
Steps to be taken in case Material is Released or Spilled Ventilate area and wash spill site			
Waste Disposal Method Mix material with combustible solvent and burn in a chemical incinerator equipped with afterburner and scrubber. Check local and state regulations. Keep tightly closed. Store in cool, dry place			
Precautions to be taken in Handling and Storing Use in chemical fume hood with proper protective lab gear.			
Other Precautions None			
Section VIII - Control Measures			
Respiratory Protection (Specify Type) MIOHSHA approved, SCBA			
Ventilation	Local Exhaust	Special	
Mechanical (General)	Required	Other	
Protective Gloves	Rubber	Eye Protection	Chem. safety goggles
Other Protective Clothing or Equipment Rubber boots			
Work/Hygiene Practices Use in chemical fume hood with proper protective lab gear.			

 Material Safety Data Sheet May be used to comply with OSHA's Hazard Communication Standard, 29 CFR 1910.1200. Standard must be consulted for specific requirements.			
IDENTITY (As Used on Label and List) InstaStain® Blue, FlashBlue™			
<small>Note: Blank space for information that is not applicable, or no information is available, this space must be marked to indicate that.</small>			
Section I Manufacturer's Name EDVOTEK, Inc. Address (Number, Street, City, State, Zip Code) 14676 Rothgeb Drive Rockville, MD 20850			
Emergency Telephone Number (301) 251-5990 Telephone Number for information (301) 251-5990 Date Prepared 03-26-09 Signature of Preparer (optional)			
Section II - Hazardous Ingredients/Identify Information Hazardous Components Specific Chemical Identity, Common Name(s) OSHA PEL ACGIH TLV Other Limits Recommended % (Optional)			
Methylene Blue Chloride No data available 3,7-Bis (Dimethylamino) Phenothiazin 5 UM Chloride No data available			
CAS # 617-334			
Section III - Physical/Chemical Characteristics			
Boiling Point	No data	Specific Gravity (H ₂ O = 1)	No data
Vapor Pressure (mm Hg.)	No data	Melting Point	No data
Vapor Density (AIR = 1)	No data	Evaporation Rate (Butyl Acetate = 1)	No data
Solubility in Water Soluble - cold			
Appearance and Odor Chemical bound to paper, no odor			
Section IV - Physical/Chemical Characteristics			
Flash Point (Method Used)	No data available	Flammable Limits	LEL No data UEL No data
Extinguishing Media Water spray, carbon dioxide, dry chemical powder, alcohol or polymer foam			
Special Fire Fighting Procedures Self contained breathing apparatus and protective clothing to prevent contact with skin and eyes			
Unusual Fire and Explosion Hazards Emits toxic fumes under fire conditions			

Section V - Reactivity Data			
Stability	Unstable Stable	Conditions to Avoid	None
Incompatibility	Strong oxidizing agents		
Hazardous Decomposition or Byproducts Toxic fumes of Carbon monoxide, Carbon dioxide, nitrogen oxides, sulfur oxides, hydrogen, chloride gas			
Hazardous Polymerization	May Occur Will Not Occur	Conditions to Avoid	None
Section VI - Health Hazard Data			
Route(s) of Entry:	Inhalation? Yes	Skin? No	Ingestion? Yes
Health Hazards (Acute and Chronic) Eyes: May cause eye irritation			
Skin: May cause skin irritation			
Eyes: May cause eye irritation			
Inhalation: Cyanosis			
NTP? IARC Monographs? OSHA Regulation? Meets criteria for proposed OSHA medical records rule PEREAC 47.30420.82			
Signs and Symptoms of Exposure No data available			
Medical Conditions Generally Aggravated by Exposure No data available			
Emergency First Aid Procedures Treat symptomatically			
Section VII - Precautions for Safe Handling and Use			
Steps to be taken in case Material is Released or Spilled Ventilate area and wash spill site			
Waste Disposal Method Mix material with a combustible solvent and burn in chemical incinerator equipped with afterburner and scrubber. Check local and state regulations. Keep tightly closed. Store in cool, dry place			
Precautions to be taken in Handling and Storing Use in chemical fume hood with proper protective lab gear.			
Other Precautions None			
Section VIII - Control Measures			
Respiratory Protection (Specify Type) MIOHSHA approved, SCBA			
Ventilation	Local Exhaust	Special	
Mechanical (General)	Required	Other	
Protective Gloves	Rubber	Eye Protection	Chem. safety goggles
Other Protective Clothing or Equipment Rubber boots			
Work/Hygiene Practices Use in chemical fume hood with proper protective lab gear.			

Mapping of Restriction Sites on Plasmid DNA

EDVOTEK Series 100 Electrophoresis Experiments:

Cat. #	Title
101	Principles and Practice of Agarose Gel Electrophoresis
102	Restriction Enzyme Cleavage Patterns of DNA
105	Mapping of Restriction Sites on Plasmid DNA
104	Size Determination of DNA Restriction Fragments
105	Mapping of Restriction Sites on Plasmid DNA
105	DNA Fingerprinting - Identification of DNA by Restriction Fragmentation Patterns
112	Analysis of <i>Eco</i> RI Cleavage Patterns of Lambda DNA
114	DNA Paternity Testing Simulation
115	Cancer Gene Detection
116	Sickle Cell Gene Detection (DNA-based)
117	Detection of Mad Cow Disease
118	Cholesterol Diagnostiics
124	DNA-based Screening for Smallpox
130	DNA Fingerprinting - Amplification of DNA for Fingerprinting



Visit our web site for information about the above experiments and other products in EDVOTEK's comprehensive offerings for biotechnology and biology education.