



The Biotechnology Education Company ®

EDVO-Kit
110
**Molecular
Weight Determination
of Proteins**

See Page 3 for storage instructions.

EXPERIMENT OBJECTIVE:

The objective of this experiment module is to determine the molecular weight of a protein using SDS horizontal gel electrophoresis. Students will develop a basic understanding of protein structure and denaturation.

EDVOTEK, Inc. • 1-800-EDVOTEK • www.edvotek.com

Molecular Weight Determination of Proteins

Table of Contents

	Page
Experiment Components	3
Experiment Requirements	3
Background Information	4
Experiment Procedures	
Experiment Overview and General Instructions	7
Preparing for Agarose Gel Electrophoresis	8
Conducting Agarose Gel Electrophoresis	13
Staining the Gel (Optional)	14
Size Determination of DNA Restriction Fragments	15
Study Questions	17
Instructor's Guidelines	
Notes to the Instructor and Pre-Lab Preparations	19
Quick Reference Tables	22
Avoiding Common Pitfalls	23
Experiment Results and Analysis	24
Study Questions and Answers	25
Material Safety Data Sheets	27

All components are intended for educational research only. They are not to be used for diagnostic or drug purposes, nor administered to or consumed by humans or animals.

THIS EXPERIMENT DOES NOT CONTAIN HUMAN DNA. None of the experiment components are derived from human sources.

EDVOTEK, The Biotechnology Education Company, and InstaStain are registered trademarks of EDVOTEK, Inc.. Ready-to-Load, Protein Agarose and LyphoProtein are trademarks of EDVOTEK, Inc.



Experiment Components

Quick Reference:

There is enough sample for 6 gels if you are using an automatic micropipet for sample delivery. Use of transfer pipets will yield fewer gels.

Reagent quantities for 6 gels are based upon the use of horizontal gel electrophoresis apparatus, Model #M12, using split trays or both halves of the standard tray.

Although the proteins in this experiment are pre-stained and can easily be visualized directly during and after electrophoresis, staining with EDVOTEK Protein InstaStain® may enhance the visibility of the bands.

READY-TO-LOAD PROTEIN SAMPLES FOR ELECTROPHORESIS

- A Pre-stained LyphoProtein™ Gel Marker (Molecular Weight Standard Protein Markers)
- B Unknown Pre-stained LyphoProtein™
- C Unknown Pre-stained LyphoProtein™
- D Unknown Pre-stained LyphoProtein™

REAGENTS & SUPPLIES

- Practice Gel Loading Solution
- Protein Agarose™ powder
- Tris-Glycine-SDS electrophoresis buffer (10x)
- 1 ml pipet
- 100 ml graduated cylinder (packaging for samples)
- Microtipped Transfer Pipets
- Semi-log graph paper template
- Protein InstaStain®

Store proteins (A-D) at -20°C.

This experiment contains ready-to-load protein samples and reagents sufficient for 6 gels (see Quick Reference).

Requirements

- Horizontal gel electrophoresis apparatus
- D.C. power supply
- Waterbath
- Recommended equipment:
 - Visualization system (white light)
 - Automatic micropipets with tips
- Pipet pump
- 250ml flasks
- Hot gloves
- Metric rulers
- Distilled or deionized water
- Large beaker
- Methanol (optional)
- Glacial acetic acid (optional)
- Small plastic tray or large weigh boat (optional)

Background Information

Proteins are a highly diversified class of biomolecules. Differences in their chemical properties, such as charge, shape, size and solubility, enables them to perform their individual, unique biological functions. These functions include enzyme catalysis, metabolic regulation, binding and transport of small molecules, gene regulation, immunological defense and cell structure. Determination of the molecular weight of a protein is of fundamental importance to its biochemical characterization. If the amino acid composition or sequence is known, the exact molecular weight can be calculated. This assumes the protein does not contain any "non-amino acid" groups (heme, zinc, covalently bonded carbohydrate, etc.) or the amount of these groups, if present, are already known. Other techniques for the determination of very accurate molecular weights include analytical ultracentrifugation and light scattering. However, these methods require large amounts of highly purified proteins and costly, sophisticated equipment. As an alternative, SDS gel electrophoresis is an easy and inexpensive method commonly used to obtain reliable estimates of protein molecular weights.

A protein can have a net negative or net positive charge, depending on its amino acid composition and specified pH environment. At certain values of pH, the molecule can be overall electrically neutral, i.e. negative and positive charges are balanced. In this case, the protein is isoelectric where the protein will not move in the electric field during electrophoresis. In the presence of an electrical field, a protein with a net charge will migrate towards the electrode of opposite charge.

Proteins exhibit many different three-dimensional shapes and folding patterns which are determined by their amino acid sequence and intracellular processing. The precise three-dimensional configuration of a protein is critical to its function. The general shapes these molecules can have are spherical, elliptical or rod-like. Proteins can consist of a single polypeptide or several polypeptides specifically associated with each other. These polypeptides can be identical, similar or completely different from one another. The number and nature of polypeptides in a protein has large effects on its mass, size and shape. Proteins that are in their normal, biologically active forms are called native.

The physical-chemical properties of proteins affect the way they migrate during gel electrophoresis. Gels used in electrophoresis (e.g. agarose, polyacrylamide) consist of microscopic pores of a defined size range that act as a molecular sieve. Only molecules with a net charge will migrate through the gel when they are in an electric field. Small molecules pass through gel pores more easily than large ones. Molecules having more of the same charge (positive or negative) than others with the same shape and size will migrate faster. Molecules of the same mass and charge can have different shapes. In this case, those with a more compact shape, like a sphere, will migrate through the gel more rapidly than those with an elongated shape, like a rod. In summary, the amount and charge (net positive or negative), the size and shape of a native protein all affect its electrophoretic migration rates. Electrophoresis of native proteins is useful in the clinical and immunological analysis of complex biological samples, such as serum.



Background Information

DENATURATION OF PROTEIN STRUCTURES

Sodium dodecylsulfate (SDS) is a detergent which consists of a hydrocarbon chain bonded to a highly negatively charged sulfate group (Figure 1). SDS binds strongly to most proteins and causes them to unfold to a random, rod-like chain. No covalent bonds are broken in this process. Therefore, the amino acid composition and sequence remains the same. Since its specific three-dimensional shape is abolished, the protein no longer possesses biological activity. Proteins that have lost their specific folding patterns and biological activity but have intact polypeptide chains, are called denatured. Proteins which contain several polypeptide chains that are associated only by non-covalent forces will be dissociated by SDS into separate, denatured polypeptide chains. Proteins can contain covalent disulfide crosslinks bonds. These bonds are formed between two cysteine amino acid residues that can be located within the same or different polypeptide chains. High concentrations of reducing agents, such as 2-mercaptoethanol, can break disulfide bonds. This allows SDS to completely dissociate and denature the protein. Figure 2 illustrates a protein containing two different sized polypeptide chains that are cross-linked by a disulfide bond. The chains are also associated by non-covalent forces. The circles represent the native structure.





Figure 1:
The chemical structure of sodium dodecylsulfate (SDS).

Figure 2:

Protein denaturation in the presence of 2-mercaptoethanol.



Background Information

In most cases, SDS binds to proteins in a constant ratio of 1.4 grams of SDS per gram of protein. On average, the number of bound SDS molecules is half the number of amino acid residues in the polypeptide. The total negative charge of SDS is much more than the intrinsic charges of the protein. SDS efficiently masks these intrinsic charges, consequently, SDS denatured proteins are net negative since the binding of the detergent is proportional to the mass of the protein. The shapes of SDS denatured proteins are the same (rod-like). The larger the molecular weight of the protein, the longer the rod-like chain. During denaturing SDS gel electrophoresis, proteins migrate through the gel towards the positive electrode at a rate that is inversely proportional to their molecular weights. In other words, the smaller the "rod-like" denatured protein, the faster it migrates. The molecular weight of the unknown is obtained by the comparison of its position after electrophoresis to the positions of standard SDS denatured proteins electrophoresed in parallel. The molecular weights of the standard proteins have been previously determined with great accuracy by some of the methods discussed previously. After the proteins have been visualized by staining the gel with dyes, their migration distance can be measured. The \log_{10} of the molecular weights of the standard proteins are plotted versus their migration distance, or RF. Taking the logarithm allows the data to be plotted as a straight line. The molecular weight of the unknown is then easily calculated from the standard curve.

For convenience, the proteins in this experiment have been partially "pre-stained". Protein samples have been covalently conjugated to a blue dye via certain amino acid side chains. This process causes the denaturation of the proteins. However, SDS is still required to produce the electrophoretic mobilities necessary for molecular weight determinations. The mobilities of the conjugates are different from those of the unconjugated forms in the presence of SDS. For this reason, the pre-stained standards cannot be reliably used to determine the molecular weight of unconjugated proteins run in parallel (visualized by staining the gel). The extent of labeling with dye can be variable, producing lot specific variances in apparent pre-stained molecular weights.

The protein standards in this experiment are a mixture having the following denatured molecular weights: 94,000; 67,000; 43,000; 30,000; 20,100. Note that these values are for the dye-conjugated proteins. The molecular weight values have been rounded off for convenience in graphical analysis.




Experiment Overview and General Instructions

EXPERIMENT OBJECTIVE:

The objective of this experiment module is to determine the molecular weight of a protein using SDS horizontal gel electrophoresis. Students will develop a basic understanding of protein structure and denaturation.

LABORATORY SAFETY

1. Gloves and goggles should be worn routinely as good laboratory practice.
2. Exercise extreme caution when working with equipment that is used in conjunction with the heating and/or melting of reagents.
3. DO NOT MOUTH PIPET REAGENTS - USE PIPET PUMPS.
4. Exercise caution when using any electrical equipment in the laboratory.
5. Always wash hands thoroughly with soap and water after handling reagents or biological materials in the laboratory.



LABORATORY NOTEBOOK RECORDINGS:

Address and record the following in your laboratory notebook or on a separate worksheet.

Before starting the Experiment:

- Write a hypothesis that reflects the experiment.
- Predict experimental outcomes.

During the Experiment:

- Record (draw) your observations, or photograph the results.



Following the Experiment:

- Formulate an explanation from the results.
- Determine what could be changed in the experiment if the experiment were repeated.
- Write a hypothesis that would reflect this change.

Preparations for Protein Agarose Gel Electrophoresis

PREPARING THE GEL BED

1. Make sure the half gel bed is clean and dry.
2. Close off the open ends of the bed by using rubber dams or tape.
 - A. Using Rubber dams:
 - Place a rubber dam on each end of the bed. Make sure the rubber dam sits firmly in contact with the sides and bottom of the bed.
 - B. Taping with labeling or masking tape:
 - With 3/4 inch wide tape, extend the tape over the sides and bottom edge of the bed.
 - Fold extended edges of the tape back onto the sides and bottom. Press contact points firmly to form a good seal.
3. Place the well forming template (comb) across the bed in the first set of notches. The comb should sit firmly and evenly across the bed of the half gel.



Caution!

This high percent agarose will boil over if not monitored. Keep a close eye on the agarose as it is being melted.

CASTING THE AGAROSE GEL

The protein agarose gel concentration required for this experiment is 3.2% weight by volume. Refer to **Table A** for guidelines.

1. Use a 400ml flask to prepare the agarose solution.

A large beaker is required to prevent the solution from boiling over during heating.

- Using a beaker, measure and add the distilled water as indicated in **Table A**. Do not add 10x buffer yet at this step. The salts in the buffer tend to cause the agarose solution to boil over.

Table A: Guidelines for Preparing Individual 3.2% Protein Agarose Gels

EDVOTEK Model #	Approximate Gel Bed Dimensions (W x L)	Amt of Protein Agarose	Volume of Distilled Water	Melt agarose in distilled water before adding 10 x buffer.	Volume of Buffer (10x)	Total Volume
M6 or M36	7 x 7 cm	0.96 gm	27 ml		3.0 ml	30 ml
M12	7 x 14 cm	1.92 gm	54 ml		6.0 ml	60 ml




Preparations for Protein Agarose Gel Electrophoresis

2. Add the required amount of Protein Agarose powder. Swirl to disperse clumps.
3. With a marking pen, indicate the level of the solution volume on the outside of the flask.
4. Heat the mixture to dissolve the agarose powder. The final solution should appear transparent.
 - A. Microwave method:
 - Cover flask with plastic wrap to minimize evaporation.
 - Heat the mixture on High for 1 minute.
 - Gently swirl the mixture and heat on high until all the agarose is dissolved. Use short bursts to avoid boiling over.
 - The mixture should become viscous and transparent with little or no particles of undissolved agarose.
 - B. Hot plate or burner method:
 - Cover the flask with foil to prevent excess evaporation.
 - Place the flask in a large beaker filled with water to avoid burning the slurry by direct exposure to heat.
 - Heat the mixture to boiling over a burner with occasional gentle swirling. Boil until all the agarose is dissolved.


DO NOT ALLOW AGAROSE TO BOIL OVER.**DO NOT LEAVE THE MIXTURE UNATTENDED DURING HEATING.**

This agarose causes rapid and vigorous boiling. Initially, the mixture will form a large head of fine bubbles that will rise quickly. When this happens, use hot gloves and remove the beaker from the heat. Let the bubbles subside, swirl briefly, and heat the beaker again. Do this several times until the mixture boils evenly.

5. Cool the protein agarose solution to 70°C with occasional gentle swirling to promote even dissipation of heat. If detectable evaporation has occurred, add distilled water to bring the solution up to the original volume as marked on the flask in step 3.



Wear gloves
and safety goggles



Useful Hint!

At high altitudes, it is necessary to use a microwave oven to reaching boiling temperatures.

After adding SDS, avoid vigorous swirling to minimize bubbling.



DO NOT POUR BOILING HOT AGAROSE INTO THE BED.

Hot agarose solution may irreversibly warp the bed.

Preparations for Protein Agarose Gel Electrophoresis

After the gel is cooled to 70°C:


6. Add the buffer concentrate (10x) as specified in Table A to the cooled agarose solution.
7. Gently swirl to mix. Try to avoid creating excess bubbles.

If using rubber dams, go to step 9. If using tape, continue with step 8.

8. Seal the interface of the bed and tape to prevent the agarose solution from leaking.
 - Use a transfer pipet to deposit a small amount of cooled agarose.
 - Wait approximately 2-3 minutes for the agarose to solidify.
9. Pour the cooled agarose solution into the bed. Make sure the bed is on a level surface.
10. Allow the gel to completely solidify. It will become firm and cool to the touch.

Remember!


The same 10x concentrated buffer is used for preparing the agarose gel buffer and the chamber buffer. SDS is added to the gel after it is melted.



Preparations for Protein Agarose Gel Electrophoresis

PREPARING THE SOLIDIFIED GEL FOR ELECTROPHORESIS

- Carefully remove the rubber dams or tape.
- Remove the comb by slowly pulling it straight up. Do this carefully and evenly to prevent tearing the sample wells. The comb may be more difficult to remove than lower percentage agarose gels generally prepared for DNA electrophoresis.
- Inspect the wells by viewing the gel from the edge nearest the wells. If some of the wells are ripped through their bottoms or sides, do not use them when loading samples.
- Place the gel in the electrophoresis chamber, properly oriented, centered and level on the platform.
- Fill the chamber of the electrophoresis apparatus with the required volume of diluted buffer as outlined in **Table B**.
- The gel is now ready for electrophoresis.



NOTE:

This 3.2% Protein Agarose gel is very sturdy and hard.

Table B: Electrophoresis (Chamber) Buffer

EDVOTEK Model #	Concentrated Buffer (10x)	Distilled Water	=	Total Volume
M6 +	30 ml	270 ml		300 ml
M12	40 ml	360 ml		400 ml
M36	100 ml	900 ml		1000 ml

Practice Gel Loading

Quick Reference:


If you are using an automatic micropipet, deliver 20 microliters to the sample well. If using transfer pipets, load the sample well until it is full.

Remember!

Use plain agar and water to make practice gels - save the agarose for the experiment.

EDVOTEK® experiments which involve electrophoresis contain practice gel loading solution. If your students are unfamiliar with loading samples in agarose gels, it is suggested that they practice sample delivery techniques before performing the electrophoresis part of an experiment. Using the EDVOTEK system, sample delivery can be performed by using either an automatic micropipet, or disposable microtipped transfer pipets.

Casting of a separate practice gel is highly recommended. One suggested activity for practice gel loading is outlined below:



1. Cast a gel with the maximum number of wells and place it under the buffer in an electrophoresis apparatus chamber. (Use plain agar and water to make the practice gels- save the agarose for the experiment.)
2. Let students practice delivering the practice gel loading solution to the sample wells.
3. If students need more practice, remove the practice gel loading solution by squirting buffer into the wells with a transfer pipet.
4. When students are finished practicing, replace the practice gel with a fresh gel and continue with the experiment.


The practice gel loading solution will become diluted in the buffer and will not interfere with the experiment.

Conducting Agarose Gel Electrophoresis

This experiment requires a 3.2% Protein Agarose gel. Make sure the electrophoresis apparatus leads reach the power source before loading samples. Do not move the apparatus after the samples are loaded because movement of the unit will cause samples to spill out of the wells.

LOADING PRE-STAINED PROTEIN SAMPLES

1. Bring a beaker of water, covered with aluminum foil, to a boil. Remove from heat.
2. Make sure resuspended tubes are tightly capped and thawed. The bottom of the tubes should be pushed through the foil and immersed in the boiling water for 3-4 minutes. The tubes should be kept suspended by the foil.
3. While the samples are still warm, load 20 μ l of each sample in tubes A - D into the wells in consecutive order.



RUNNING THE GEL

1. After the samples are loaded, carefully snap the cover down onto the electrode terminals.

Make sure that the negative and positive indicators on the cover and apparatus chamber are properly oriented.
2. Insert the plug of the black wire into the black input of the power source (negative input). Insert the plug of the red wire into the red input of the power source (positive input).
3. Set the power source at the required voltage and run the electrophoresis for the length of time as determined by your instructor. See **Table C** for time and voltage guidelines.

When current is flowing properly, you should see bubbles forming on the electrodes.

4. After the electrophoresis is completed, turn off the power, unplug the power source, disconnect the leads and remove the cover.
5. Transfer the gel to a white light box to enhance visualization of the pre-stained protein bands.
6. If the measurements cannot be done shortly after the electrophoresis is completed, you may wish to stain and/or preserve the gel by performing the steps which follow.

Reminders:

During electrophoresis, the protein samples migrate through the agarose gel towards the positive electrode. Before loading the samples, make sure the gel is properly oriented in the apparatus chamber.





Table C: Time and Voltage

Volts	Recommended Time	
	Minimum	Optimal
125	30 min	45 min
70	40 min	1.5 hrs
50	60 min	2.0 hrs

Staining the Gel (Optional)

Although the proteins in this experiment are pre-stained and can easily be visualized directly during and after electrophoresis, staining with EDVOTEK® Protein InstaStain® may enhance the visibility of the bands.

STAINING WITH PROTEIN INSTASTAIN® IN ONE EASY STEP

Protein agarose gels can be stained with Protein InstaStain® cards in one easy step.

1. After electrophoresis, submerge the gel in a small tray with 100 ml of fixative solution. (Use enough solution to cover the gel.)
2. Gently float a sheet of Protein InstaStain® with the stain side (blue) in the liquid. Cover the gel to prevent evaporation.
3. Gently agitate on a rocking platform for 1-3 hours or overnight.
4. After staining, Protein bands will appear as dark blue bands against a light background and will be ready for photography.
NO DESTAINING IS REQUIRED.
5. If the gel is too dark, destain in several changes of fresh destain solution until the appearance and contrast of the protein bands against the background improves.

Storing the Gel

For long term storage, the gel should be stored in a mixture of 50 ml of distilled water containing 6 ml of acetic acid and 3 ml of glycerol overnight.

**Fixative/Destaining
Solution for each gel
(100ml)**

50 ml	Methanol
10 ml	Glacial Acetic Acid
40 ml	Distilled Water



Size Determination of DNA Restriction Fragments

The first step for determining the molecular weight of the "unknown" protein polypeptides is to determine the migration distances of the fragments generated after electrophoresis. The molecular weight of the "unknowns" will be extrapolated by graphing their migration distances relative to Sample A, Prestained Lyophilized Protein Gel Markers (Standard markers), for which the molecular weight of each fragment is known.


1. Measure and record the distance traveled in the agarose gel by each Standard marker. In each case, measure from the lower edge of the sample well to the lower end of each band. Record the distance traveled in centimeters (to the nearest millimeter).
2. Label the semi-log graph paper:
 - A. Label the non-logarithmic horizontal x-axis "Migration Distance" in centimeters at equal intervals.
 - B. Label the logarithmic vertical y-axis "Molecular Weight".


Choose your scales so that the data points are well spread out. Assume the first cycle on the y-axis represents 1,000 to 10,000, the second cycle represents 10,000 to 100,000, and the third cycle 100,000 to 1,000,000.

3. For each Standard marker, plot the measured migration distance on the x-axis versus its molecular weight on the y-axis.
4. Draw the best average straight line through all the points.

The line should have approximately equal numbers of points scattered on each side of the line. Some points may be right on the line (see example at left).
5. Measure the migration distance of each of the "unknown" proteins.
6. Using the graph of the Standard markers, determine the molecular weight of each "unknown".
7. Find the migration distance of the unknown protein on the x-axis. Then draw a vertical line from that point until the standard graph line is intersected.
8. From the point of intersection, draw a second line horizontally to the y-axis and determine the apparent molecular weight of the protein polypeptide.

The assignment of molecular weight for proteins separated by agarose gel electrophoresis has a $\pm 10\%$ margin of error.







Material Safety Data Sheets

Full-size (8.5 x 11") pdf copy of MSDS is available at www.edvotek.com or by request.


110
Experiment

 <p>Material Safety Data Sheet May be used to comply with OSHA's Hazard Communication Standard. 29 CFR 1910.1200 Standard must be consulted for specific requirements.</p>			
IDENTITY (As Used on Label and List) Sodium Dodecyl Sulfate (SDS)		Note: Blank spaces are not permitted. If any item is not applicable, or no information is available, the space must be marked to indicate that.	
Section I Manufacturer's Name EDVOTEK, Inc. Address (Number, Street, City, State, Zip Code) 14676 Rothgeb Drive Rockville, MD 20850		Emergency Telephone Number (301) 251-5990 Telephone Number for information (301) 251-5990 Date Prepared 08/04/09 Signature of Preparer (optional)	
Section II - Hazardous Ingredients/Identify Information Hazardous Components (Specific Chemical Identity; Common Name(s)) OSHA PEL ACGIH TLV Other Limits Recommended % (Optional) Lauryl Sulfate, Sodium No data No data No data No data C12H26O4S CAS# 151-21-3			
Section III - Physical/Chemical Characteristics			
Boiling Point	No data	Specific Gravity (H ₂ O = 1)	No data
Vapor Pressure (mm Hg.)	No data	Melting Point	No data
Vapor Density (AIR = 1)	No data	Evaporation Rate (Butyl Acetate = 1)	No data
Solubility in Water		Soluble	
Appearance and Odor		Clear liquid, no odor	
Section IV - Physical/Chemical Characteristics			
Flash Point (Method Used)	No data	Flammable Limits	LEL No data UEL No data
Extinguishing Media Water spray, carbon dioxide, dry chemical powder, alcohol or polymer foam			
Special Fire Fighting Procedures Wear SCBA and protective clothing to prevent contact with skin & eyes.			
Unusual Fire and Explosion Hazards May emit toxic fumes.			


Section V - Reactivity Data			
Stability	Unstable	Conditions to Avoid	
	Stable	X	None
Incompatibility Strong oxidizing agents			
Hazardous Decomposition or Byproducts Carbon monoxide, carbon dioxide, sulfur oxides			
Hazardous Polymerization	May Occur	Conditions to Avoid	
	Will Not Occur	X	None
Section VI - Health Hazard Data			
Route(s) of Entry:		Inhalation?	Skin? Ingestion?
		Yes	Yes Yes
Health Hazards (Acute and Chronic) May cause irritation to eyes, ears and nose.			
Carcinogenicity:		NTP?	IARC Monographs? OSHA Regulation?
		No data	No data
Signs and Symptoms of Exposure Respiratory tract: burning sensation, coughing, wheezing, laryngitis, shortness of breath, & headache			
Medical Conditions Generally Aggravated by Exposure No data			
Emergency First Aid Procedures Flush skin/eyes with large amounts of water. If inhaled, remove to fresh air.			
Section VII - Precautions for Safe Handling and Use			
Steps to be Taken in case Material is Released for Spilled Evacuate area. Wear SCBA, rubber boots and rubber gloves. Mop up with absorbent material and burn in chemical incinerator equipped with an afterburner and scrubber.			
Waste Disposal Method Observe all federal, state, and local laws.			
Precautions to be Taken in Handling and Storing Wear protective gear. Avoid contact/inhalation.			
Other Precautions Strong sensitizers			
Section VIII - Control Measures			
Respiratory Protection (Specify Type) NIOSH/MSHA approved respirator.			
Ventilation	Local Exhaust	No	Special Chem. fume hood
	Mechanical (General)	No	Other None
Protective Gloves	rubber		Eye Protection Splash proof goggles
Other Protective Clothing or Equipment Rubber boots			
Work/Hygienic Practices Avoid prolonged or repeated exposure.			

 <p>Material Safety Data Sheet May be used to comply with OSHA's Hazard Communication Standard. 29 CFR 1910.1200 Standard must be consulted for specific requirements.</p>			
IDENTITY (As Used on Label and List) Protein Agarose		Note: Blank spaces are not permitted. If any item is not applicable, or no information is available, the space must be marked to indicate that.	
Section I Manufacturer's Name EDVOTEK, Inc. Address (Number, Street, City, State, Zip Code) 14676 Rothgeb Drive Rockville, MD 20850		Emergency Telephone Number (301) 251-5990 Telephone Number for information (301) 251-5990 Date Prepared 08/04/09 Signature of Preparer (optional)	
Section II - Hazardous Ingredients/Identify Information Hazardous Components (Specific Chemical Identity; Common Name(s)) OSHA PEL ACGIH TLV Other Limits Recommended % (Optional) This product contains no hazardous materials as defined by the OSHA Hazard Communication Standard.			
Section III - Physical/Chemical Characteristics			
Boiling Point	65°C	Specific Gravity (H ₂ O = 1)	No data
Vapor Pressure (mm Hg.)	No data	Melting Point	No data
Vapor Density (AIR = 1)	No data	Evaporation Rate (Butyl Acetate = 1)	No data
Solubility in Water		Souluble-hot, insoluble-cold	
Appearance and Odor		White powder, no odor	
Section IV - Physical/Chemical Characteristics			
Flash Point (Method Used)	Not determined	Flammable Limits	LEL No data UEL No data
Extinguishing Media Water spray, dry chemical, carbon dioxide, halon or standard foam			
Special Fire Fighting Procedures Negligible fire hazard when exposed to heat or flame			
Unusual Fire and Explosion Hazards None			

Section V - Reactivity Data			
Stability	Unstable	Conditions to Avoid	
	Stable	X	None
Incompatibility None			
Hazardous Decomposition or Byproducts None hazardous			
Hazardous Polymerization	May Occur	Conditions to Avoid	
	Will Not Occur	X	None
Section VI - Health Hazard Data			
Route(s) of Entry:		Inhalation?	Skin? Ingestion?
		Yes	Yes Yes
Health Hazards (Acute and Chronic) Inhalation: No data available Ingestion: acute large amts cause diarrhea, flatulence & rarelyfecal impaction Skin: no data			
Carcinogenicity:		NTP?	IARC Monographs? OSHA Regulation?
		None, on FDA list of direct food substances & generally recognized as safe	
Signs and Symptoms of Exposure No data available			
Medical Conditions Generally Aggravated by Exposure No data available except as listed above			
Emergency First Aid Procedures Treat symptomatically and supportively			
Section VII - Precautions for Safe Handling and Use			
Steps to be Taken in case Material is Released for Spilled Sweep up and place in suitable containers for reclamation or disposal			
Waste Disposal Method Normal solid waste disposal			
Precautions to be Taken in Handling and Storing None			
Other Precautions None			
Section VIII - Control Measures			
Respiratory Protection (Specify Type) Chem. cartridge respirator with full facepiece, organic vapor cart.			
Ventilation	Local Exhaust	Special	
	Mechanical (General)	Gen. dilution vent.	Other
Protective Gloves	Yes		Eye Protection Splash proof goggles
Other Protective Clothing or Equipment Impervious clothing to prevent skin contact			
Work/Hygienic Practices None			

 <p align="center">Material Safety Data Sheet May be used to comply with OSHA's Hazard Communication Standard. 29 CFR 1910.1200. Standard must be consulted for specific requirements.</p>				
IDENTITY (As Used on Label and List) Tris-Glycine-SDS Buffer		Note: Blank spaces are not permitted. If any item is not applicable, or no information is available, the space must be marked to indicate that.		
Section I				
Manufacturer's Name EDVOTEK		Emergency Telephone Number 301-251-5990		
Address (Number, Street, City, State, and ZIP Code) 14676 Rothgeb Drive Rockville, MD 20850		Telephone Number for information 301-251-5990		
		Date Prepared 08/04/09		
		Signature of Preparer (optional)		
Section II - Hazardous Ingredients/Identify Information				
Hazardous Components (Specific Chemical Identity; Common Name(s)) OSHA PEL ACGIH TLV Other Limits Recommended % Optional				
Lyl Sulfate, Sodium No data No data No data No data				
Section III - Physical/Chemical Characteristics				
Boiling Point	No data	Specific Gravity (H20 = 1)	No data	
Vapor Pressure (mm Hg.)	No data	Melting Point	No data	
Vapor Density (AIR =1)	No data	Evaporation Rate (Butyl Acetate =1)	No data	
Solubility in Water Soluble				
Appearance and Odor Clear liquid, no odor				
Section IV - Fire and Explosion Hazard Data N.D = No data				
Flash Point (Method Used)	No data	Flammable Limits	LEL N.D	UEL N.D
Extinguishing Media Water spray, Carbon dioxide, dry chemical powder, alcohol or polymer foam				
Special Fire Fighting Procedures Wear SCBA and protective clothing to prevent contact w/ skin and eyes.				
Unusual Fire and Explosion Hazards May emit toxic fumes.				

Section V - Reactivity Data				
Stability	Unstable		Conditions to Avoid	
	Stability	X	None	
Incompatibility (Materials to avoid) Strong oxidizing agents				
Hazardous Decomposition or Byproducts Carbon monoxide, carbon dioxide, sulfur oxides.				
Hazardous Polymerization	May Occur		Conditions to Avoid	
	Will Not Occur	X	None	
Section VI - Health Hazard Data				
Route(s) of Entry: Inhalation? Yes Skin? Yes Ingestion? Yes				
Health Hazards (Acute and Chronic) Irritating to eyes, mucous membranes and upper respiratory tract. Chronic exposure may cause lung damage or pulmonary sensitization resulting in hyperactive airway dysfunction.				
Carcinogenicity: NTP? No data IARC Monographs? No data OSHA Regulation? No data				
Signs and Symptoms of Exposure				
Respiratory tract: burning sensation, coughing, wheezing, laryngitis, shortness of breath, headache, nausea				
Medical Conditions Generally Aggravated by Exposure No data				
Emergency and First Aid Procedures				
Flush skin/eyes w/ large amounts of water. Inhaled: remove to fresh air.				
Ingestion: Give large amounts of water or milk. Do not induce vomiting.				
Section VII - Precautions for Safe Handling and Use				
Steps to be Taken in case Material is Released or spilled Evacuate area. Wear SCBA, rubber boots and rubber gloves.				
Mop up w/ absorbent material and burn in chemical incinerator equipped w/ afterburner & scrubber.				
Waste Disposal Method Observe all federal, state, and local laws.				
Precautions to be Taken in Handling and Storing Wear protective gear. Avoid contact/inhalation.				
Other Precautions Strong sensitizer				
Section VIII - Control Measures				
Respiratory Protection (Specify Type) NIOSH/MSHA approved respirator				
Ventilation	Local Exhaust	No	Special	Chem. fume hood
	Mechanical (General)	No	Other	None
Protective Gloves	Rubber		Eye Protection	Splash proof goggles
Other Protective Clothing or Equipment Rubber boots				
Work/Hygenic Practices Avoid prolonged or repeated exposure.				

 <p align="center">Material Safety Data Sheet May be used to comply with OSHA's Hazard Communication Standard. 29 CFR 1910.1200. Standard must be consulted for specific requirements.</p>				
IDENTITY (As Used on Label and List) Protein InstaStain		Note: Blank spaces are not permitted. If any item is not applicable, or no information is available, the space must be marked to indicate that.		
Section I				
Manufacturer's Name EDVOTEK, Inc.		Emergency Telephone Number (301) 251-5990		
Address (Number, Street, City, State, Zip Code) 14676 Rothgeb Drive Rockville, MD 20850		Telephone Number for information (301) 251-5990		
		Date Prepared 08/04/09		
		Signature of Preparer (optional)		
Section II - Hazardous Ingredients/Identify Information				
Hazardous Components (Specific Chemical Identity; Common Name(s)) OSHA PEL ACGIH TLV Other Limits Recommended % (Optional)				
Methanol (Methyl Alcohol) 200ppm 200ppm No data 90%-100% CH3OH				
Section III - Physical/Chemical Characteristics				
Boiling Point	65°C	Specific Gravity (H ₂ O = 1)	.79	
Vapor Pressure (mm Hg.)	96mmHg	Melting Point	N/A	
Vapor Density (AIR = 1)	1.11	Evaporation Rate (Butyl Acetate = 1)	4.6	
Solubility in Water Complete (100%)				
Appearance and Odor chemical bound to paper, no odor				
Section IV - Physical/Chemical Characteristics				
Flash Point (Method Used)	(closed cup) 12°C	Flammable Limits	LEL 6.0%	UEL 36%
Extinguishing Media Use alcohol foam, dry chemical or carbon dioxide. (Water may be ineffective)				
Special Fire Fighting Procedures Wear SCBA with full facepiece operated in positive pressure mode. Move containers from firearea				
Unusual Fire and Explosion Hazards Vapors may flow along surfaces to distant ignition sources.				


Section V - Reactivity Data				
Stability	Unstable		Conditions to Avoid	
	Stable	X	None	
Incompatibility Strong oxidizing agents				
Hazardous Decomposition or Byproducts Carbon monoxide, Carbon dioxide, Sulfur oxides				
Hazardous Polymerization	May Occur		Conditions to Avoid	
	Will Not Occur	X	None	
Section VI - Health Hazard Data				
Route(s) of Entry: Inhalation? Yes Skin? Yes Ingestion? Yes				
Health Hazards (Acute and Chronic) Irritating to eyes, skin, mucous membranes and upper respiratory tract. Chronic exposure may cause lung damage or pulmonary sensitization				
Carcinogenicity: NTP? No data IARC Monographs? No data OSHA Regulation? No data				
Signs and Symptoms of Exposure				
Respiratory tract: burning sensation. Coughing, wheezing, laryngitis, shortness of breath, headache				
Medical Conditions Generally Aggravated by Exposure No data				
Emergency First Aid Procedures				
Flush skin/eyes w/ large amounts of water. If inhaled, remove to fresh air. Ingestion: give large amounts of water or milk. Do not induce vomiting.				
Section VII - Precautions for Safe Handling and Use				
Steps to be Taken in case Material is Released for Spilled Evacuate area. Wear SCBA, rubber boots and rubber gloves. Mop up w/ absorbent material and burn in chemical incinerator equipped w/ an afterburner and scrubber.				
Waste Disposal Method Observe all federal, state, and local laws.				
Precautions to be Taken in Handling and Storing Wear protective gear. Avoid contact/inhalation.				
Other Precautions Strong sensitizer				
Section VIII - Control Measures				
Respiratory Protection (Specify Type) NIOSH/MSHA approved respirator				
Ventilation	Local Exhaust	No	Special	Chem fume hood
	Mechanical (General)	No	Other	None
Protective Gloves	Rubber		Eye Protection	Splash-proof goggles
Other Protective Clothing or Equipment Rubber boots				


Material Safety Data Sheets

Full-size (8.5 x 11") pdf copy of MSDS is available at www.edvotek.com or by request.

110

Experiment

 <p>Material Safety Data Sheet May be used to comply with OSHA's Hazard Communication Standard. 29 CFR 1910.1200 Standard must be consulted for specific requirements.</p>			
IDENTITY (As Used on Label and List) Protein Samples		Note: Blank spaces are not permitted. If any item is not applicable, or no information is available, the space must be marked to indicate that.	
Section I		Emergency Telephone Number (301) 251-5990	
Manufacturer's Name EDVOTEK, Inc.		Telephone Number for information (301) 251-5990	
Address (Number, Street, City, State, Zip Code) 14676 Rothgeb Drive Rockville, MD 20850		Date Prepared 08/04/09	
		Signature of Preparer (optional)	
Section II - Hazardous Ingredients/Identify Information			
Hazardous Components (Specific Chemical Identity: Common Name(s)) OSHA PEL ACGIH TLV Other Limits Recommended % (Optional)			
Lauryl Sulfate, Sodium No data No data No data No data			
C ₁₂ H ₂₆ O ₄ S			
CAS # 151-21-3			
Section III - Physical/Chemical Characteristics			
Boiling Point	No data	Specific Gravity (H ₂ O = 1)	No data
Vapor Pressure (mm Hg.)	No data	Melting Point	No data
Vapor Density (AIR = 1)	No data	Evaporation Rate (Butyl Acetate = 1)	No data
Solubility in Water Soluble			
Appearance and Odor Clear liquid, characteristic, disagreeable odor			
Section IV - Physical/Chemical Characteristics			
Flash Point (Method Used)	No data	Flammable Limits	LEL No data UEL No data
Extinguishing Media Water spray, carbon dioxide, dry chemical powder, alcohol or polymer foam			
Special Fire Fighting Procedures Wear SCBA and protective clothing to prevent contact with skin.			
Unusual Fire and Explosion Hazards May emit toxic fumes			

 <p>Material Safety Data Sheet May be used to comply with OSHA's Hazard Communication Standard. 29 CFR 1910.1200 Standard must be consulted for specific requirements.</p>			
IDENTITY (As Used on Label and List) Practice Gel Loading Solution		Note: Blank spaces are not permitted. If any item is not applicable, or no information is available, the space must be marked to indicate that.	
Section I		Emergency Telephone Number (301) 251-5990	
Manufacturer's Name EDVOTEK, Inc.		Telephone Number for information (301) 251-5990	
Address (Number, Street, City, State, Zip Code) 14676 Rothgeb Drive Rockville, MD 20850		Date Prepared 08/04/09	
		Signature of Preparer (optional)	
Section II - Hazardous Ingredients/Identify Information			
Hazardous Components (Specific Chemical Identity: Common Name(s)) OSHA PEL ACGIH TLV Other Limits Recommended % (Optional)			
This product contains no hazardous materials as defined by the OSHA Hazard Communication Standard.			
Section III - Physical/Chemical Characteristics			
Boiling Point	No data	Specific Gravity (H ₂ O = 1)	No data
Vapor Pressure (mm Hg.)	No data	Melting Point	No data
Vapor Density (AIR = 1)	No data	Evaporation Rate (Butyl Acetate = 1)	No data
Solubility in Water Soluble			
Appearance and Odor Blue liquid, no odor			
Section IV - Physical/Chemical Characteristics			
Flash Point (Method Used)	No data	Flammable Limits	LEL No data UEL No data
Extinguishing Media Dry chemical, carbon dioxide, water spray or foam			
Special Fire Fighting Procedures Use agents suitable for type of surrounding fire. Keep upwind, avoid breathing hazardous sulfur oxides and bromides. Wear SCBA.			
Unusual Fire and Explosion Hazards Unknown			

Section V - Reactivity Data			
Stability	Unstable	Conditions to Avoid	
	Stable	X	None
Incompatibility Strong oxidizing agents			
Hazardous Decomposition or Byproducts Carbon monoxide, carbon dioxide, sulfur oxides			
Hazardous Polymerization	May Occur	Conditions to Avoid	
	Will Not Occur	X	None
Section VI - Health Hazard Data			
Route(s) of Entry:	Inhalation?	Skin?	Ingestion?
	Yes	Yes	Yes
Health Hazards (Acute and Chronic) May cause irritation to eyes, ears and nose.			
Carcinogenicity:	NTP?	IARC Monographs?	OSHA Regulation?
		No data	
Signs and Symptoms of Exposure Respiratory tract: burning sensation, coughing, wheezing, laryngitis, shortness of breath, & headache			
Medical Conditions Generally Aggravated by Exposure No data			
Emergency First Aid Procedures Flush skin/eyes with large amounts of water. If inhaled, remove to fresh air. Ingestion: give large amounts of water or milk. Do not induce vomiting.			
Section VII - Precautions for Safe Handling and Use			
Steps to be Taken in case Material is Released for Spilled Evacuate area. Wear SCBA, rubber boots and rubber gloves. Mop up with absorptive material and burn in chemical incinerator equipped with afterburner and scrubber.			
Waste Disposal Method Observe all federal, state, and local laws.			
Precautions to be Taken in Handling and Storing Wear protective gear. Avoid contact/inhalation.			
Other Precautions Strong sensitizer			
Section VIII - Control Measures			
Respiratory Protection (Specify Type) NIOSH/MSHA approved respirator			
Ventilation	Local Exhaust	No	Special Chem. fume hood
	Mechanical (General)	No	Other None
Protective Gloves	Rubber	Eye Protection	Splash proof goggles
Other Protective Clothing or Equipment Rubber boots			
Work/Hygienic Practices Avoid prolonged or repeated exposure			

Section V - Reactivity Data			
Stability	Unstable	Conditions to Avoid	
	Stable	X	None
Incompatibility None			
Hazardous Decomposition or Byproducts Sulfur oxides, and bromides			
Hazardous Polymerization	May Occur	Conditions to Avoid	
	Will Not Occur	X	None
Section VI - Health Hazard Data			
Route(s) of Entry:	Inhalation?	Skin?	Ingestion?
	Yes	Yes	Yes
Health Hazards (Acute and Chronic) Acute eye contact: May cause irritation. No data available for other routes.			
Carcinogenicity:	NTP?	IARC Monographs?	OSHA Regulation?
Signs and Symptoms of Exposure May cause skin or eye irritation			
Medical Conditions Generally Aggravated by Exposure None reported			
Emergency First Aid Procedures Treat symptomatically and supportively. Rinse contacted area with copious amounts of water.			
Section VII - Precautions for Safe Handling and Use			
Steps to be Taken in case Material is Released for Spilled Wear eye and skin protection and mop spill area. Rinse with water.			
Waste Disposal Method Observe all federal, state, and local regulations.			
Precautions to be Taken in Handling and Storing Avoid eye and skin contact.			
Other Precautions None			
Section VIII - Control Measures			
Respiratory Protection (Specify Type)			
Ventilation	Local Exhaust	Yes	Special None
	Mechanical (General)	Yes	Other None
Protective Gloves	Yes	Eye Protection	Splash proof goggles
Other Protective Clothing or Equipment None required			
Work/Hygienic Practices Avoid eye and skin contact			